

CONTINUOUS VIBRATING SCREEN



Non contractual photo

SERVICE : POWER SUPPLY: 220 V, SINGLE PHASE, 50 HZ

REFERENCE: MP412

Before compression, a homogeneous particle size cut should be selected. With the MP412 vibrating screen, you will do this operation without difficulty.

- Study of the unit calibration operation
- · Plotting a grain size distribution curve

Technical specifications:

- · Electromagnetic motor drive,
- · Diameter of sieves, 200mm,
- 2 sieves of cut of 63mm and 1.5mm,
- Four sieves 125, 250, 530, 850 ?m,
- · Adjustment of time and amplitude,
- · Loading container with lid,
- Three drain / recovery chutes,
- · Model to ask on bench.

OPTIONS:

Other sieves on request